INcwb072\_ES\_PT

712-047\_712-047V01

ES

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## Rollator Manual

IMPORTANT, RETAIN FOR FUTURE REFERENCE: READ CAREFULLY

# ASSEMBLY & INSTRUCTION MANUAL

# SAFETY AND WARNING INSTRUCTIONS

DO NOT exceed 300 lb. weight capacity.

DO NOT use as a wheelchair or to transport someone

DO NOT self-propel or scoot around while seated.

Care should be taken to ensure that all moving parts are in good working order before operation.

Check handle push buttons to make sure they are secure. Always contact your physician or therapist before using to determine proper adjustment and usage.

Brakes must be in locked position before using the seat. Do not attempt to push yourself or ambulate in any way.

Both feet should rest flat on the ground when sitting.

DO NOT push back on the backrest or extend your body over the backrest. It is not designed to support the entire weight of the user. Backrest must always be attached when rollator is in use.

DO NOT use on stairs, escalators, moving sidewalks, or any uneven surfaces.

DO NOT attempt to reach for items while seated on the rollator. This motion may cause a change in weight distribution and cause the rollator to tip or become unstable.

The rollator bag/basket has a weight capacity of 10 lb. Items should not hang over the top or front of the basket.

Make sure rollator frame is fully opened and in a locked position before use.

DO NOT sit on the rollator on an incline.

DO NOT support your full body weight on the handles when attempting to sit down.

# PRODUCT FEATURES

- |                          |                         |
|--------------------------|-------------------------|
| ① Handgrip               | ⑧ Brake adjustment knob |
| ② Parking/hand brake     | ⑨ Front wheel           |
| ③ Armrest pad            | ⑩ Seat                  |
| ④ Height adjustment tube | ⑪ Seat rail             |
| ⑤ Height adjustment lock | ⑫ Backrest              |
| ⑥ Fold Adjustment Lock   | ⑬ Personal item bag     |
| ⑦ Rear wheel             | ⑭ Cane Holder           |

# ASSEMBLY INSTRUCTIONS


Carefully remove the accessories from the box and remove all protective-packaging.

- **Unfold the rollator base**

1. Unfold and Stand the rollator Base on all four wheels.
2. Hold the handle and push it upward to make the folding device buckle with the main frame.

3. Push down and out on the seat support tubes until they are secured in the frame cradles and you should hear the seat click when properly engaged.

# ASSEMBLY INSTRUCTIONS

## ● Handle and forearm adjustment

Press the orange release button and push the armrest tube downward to your desired height. There should be an audible "click" to signal that the orange release button is fully engaged. Repeat the steps for the other arm rest. The armrests should be set at a height that is high enough that you feel support under your forearms when standing upright, but not so high as it may cause your shoulders to "hunch."

Set the armrests at the same height and put equal weight on each forearm while using.

## CAUTION

Prior to each use, please ensure the following:

1. Both adjustment buttons are engaged with the handles in a locked position.
2. Both handles are set at the same adjustment.

Warning: Failure to follow safety instruction may cause serious injury or death.

## ● Attaching The Cane Holder

1. Take out the cane holder ( if available) from the personal bag.
2. Snap the cane holder ( if available) all the way until it's secured.

## ATTACHING / REMOVING THE PERSONAL ITEM BAG

The personal item bag has plastic clips which go over the seat rail extension on the front end of the Rollator (located at front end of the seat) , slide the clips over the right and left seat rail extensions.

## ATTACHING/ REMOVING THE BACKREST

With the logo side of the backrest upright and facing the front of the Rollator, slide the T-slot on the end of the backrest onto the respective mating T-ridge on the front of the frame under the armrests. Repeat for the other end of the backrest. Check to confirm that the T-slots and T-ridges are fully engaged on both ends of the backrest.

Backrest may be left on the Rollator or removed when not in use.

# USING BRAKES /PARKING BRAKES

To apply the brakes while walking with the Rollator, pull both brake handles evenly toward you.

To set the parking brakes, while using one hand to hold the handgrip use the other hand to push the brake handle forward until it clicks in locked position.

## ADJUSTING BRAKES

Basic Adjustments:

To tighten the brakes if too loose, locate the orange screw knob under the frame near the rear wheel . Press the knob downwards and turn to the left (counterclockwise) a half rotation at a time until desired tightness is reached.

To loosen the brakes if too tight, press the orange screw knob downwards and turn to the right (clockwise) until desired tightness is reached. Always ensure the parking brake is disengaged when adjusting brakes.

## **How to Properly Fit Your Rollator**

Rollators are not One Size Fits All, and it is very important to choose a rollator that fits your height. Two things to consider are 1(handle height for proper walking posture and 2) seat height - which should allow you to sit safely and comfortably on the rollator with your back resting on the backrest, your feet flat on the ground and your legs at a 90 degree angle.

## **How to sit with a rollator**

1. Lock the brakes by pushing down on the brake handles until it locks in the downward position. Brakes must be on the lock position for the sitting.
2. Test the rollator before sitting down to make sure that it will not roll.
3. Use the handles to stabilize your body as you turn around.
4. Before sitting, touch the back of your legs to the edge of the seat to ensure proper positioning of the rollator.
5. Holding the handles, slowly lower your body onto the seat.

DO NOT support your full body weight on the handles when attempting to sit down. The rollator is intended for stationary seating only. If someone sits then pushes with their feet, the walker could easily tip and cause a fall. The walker should not be used as a wheel- chair and it should never be pushed or rolled when someone is sitting on it.

## **How to stand with a rollator**

1. Make sure the brakes are engaged in the downward position
2. Sit as comfortably on the edge of your seat
3. Keep your toes in line with the edge of your seat
4. Place both hands on the chair or one hand on the chair and one hand on the rollator
5. Never put all of your weight on the rollator, as it could make the rollator move
6. Lean forward until some of your body weight is on your feet
7. Use your legs to stand, your arms should only be lifting the weight your legs cannot
8. Do not move forward until you are certain about your balance
9. Unlock brakes and walk.

## How to walk with a rollator

1. Move between the arm pads, while standing in an upright position with good posture.
2. Hold the brake handles with your hands while your forearms rest on the platform.
3. This position will allow you to quickly stop the rollator or control your speed if necessary by squeezing the handbrake.
4. Roll your rollator slowly forward so that the rear wheels are a few inches ahead of your body.
5. Step forward and place one foot in line with the rear wheels.
6. Maintaining a steady roll forward, place your other foot in line with the rear wheels.
7. Determine a pace that is comfortable and allows you to place each step in line with the rear wheels.

## Maintenance

Check your rollator periodically to make sure that the brakes are working correctly and that all of the nuts and bolts are secure. DO NOT use the rollator if the brakes are not working correctly.

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