

Technical Parameters:

Rated input voltage: AC 100-220V 50-60HZ
DC 12V 800MA

Rated power: < 20W

Tai Chi Massage Mattress

User Manual



Working Principle

Thanks to your purchase our series micro-computer massage cushion, this product using the latest high-tech micro-computer technology, combining ergonomics with principles of traditional Chinese medicine massage. It sets a variety of professional masseurs massage (such as scraping, kneading, acupressure, tapping, rolling, massage, vibration), benefitting for human head, shoulders, waist, thighs, and other parts, subject to a variety of massage and relax the body's bones, reducing the body's tension and pressure, and enhance the capacity of tissue metabolism, promote blood circulation, improve the capacity of joint activities, back massage can improve the body's immune, especially for middle-age and elderly healthcare to balance Yin and Yang, promote blood circulation and improve organ function, resulting in improved body strength. Let your massage enjoy.

Operation

- ★ **Power:** Press "ON/OFF" key, begin or stop the massage function, circulation press this button, can choose working time: "10 mins", "20 mins" or "30 mins".
- ★ **Mode:** Press "Mode" key, choose 8 different working program, shows "F1-F8" on the display window.
- ★ **Strong/Weak:** Press "Strong" / "Weak" key to increase/weaken the massage strength, shows "S1-S5" on the display window.
- ★ **Fast/Slow:** Press "Fast" / "Slow" key to control speed, shows "H1-H5" on the display window.
- ★ **Infrared heat:** Press "Infrared Heat" key to open the infrared heat function.



Hand Control Schematic Diagram

Points Of Attention

- ★ After use, or before cleaning must unplug the power plug.
- ★ The machine must not be used in damp environment and in the bathroom, please use after stored in a dry and no direct sunlight place.
- ★ Do not use in sleep, to prevent irritation long time on a part of the health effects.
- ★ The use of massage machine can not be after half an hour, or because the body physique is different and the adverse effect.
- ★ In the process of using excessive stimulation and cause discomfort, should immediately stop using.
- ★ When a child or body not free when using this product, must have the right side guidance counseling.
- ★ During pregnancy, the spine bone pain, acute disease, heart disease and other body had abnormal condition, should be used under the guidance of a doctor.
- ★ Please 500CC drink boiled water or hot tea after the massage, to promote the excretion of metabolic substances.
- ★ The machine power supply cord once damaged, must consult the factory power supply special household cord exchange, or with the manufacturer, must not be to exchange unqualified cord.
- ★ This machine is in use such as surface heating, heat is not sensitive to the people when in use must pay attention to.

Analysis And Treatment Of Simple Fault

The common fault phenomena	The possible reason	Processing method
The whole machine is made completely	<ul style="list-style-type: none">★ If power failure;★ Whether the adapter into a power outlet;★ The power switch is open;★ Adapter DC plug and whether the car mat DC socket on the plug;★ The cigarette lighter is inserted into the cigarette lighter socket.	<ul style="list-style-type: none">★ After the call to use;★ Correct securely adapter into an AC power supply socket;★ Turn on the power switch;★ Strong DC plug and socket on the plug DC;★ The correct cigarette lighter into car cigarette lighter socket.
The noise	If the right place	The position is adjusted.
Too weak	<ul style="list-style-type: none">★ Network voltage is too low;★ The cigarette lighter and whether the car cigarette lighter socket on the plug firmly.	<ul style="list-style-type: none">★ Use the regulator;★ Pull out the cigarette lighter plug and then wipe the car cigarette lighter socket on the plug.